

# miniRaman Spectrometer

The world's smallest  
Raman Spectrometer



- Premium performance for quantitative analysis
- Excellent for material identification
- Accuracy - high sensitivity
- Extended spectral range via dual laser
- Affordable

# Software

## Miraspec for PC (Windows 10, 11) and smartphone (Android)

System controlled by smartphone or PC via Bluetooth or USB-C cable.

|                                   |   |  |
|-----------------------------------|---|--|
| <b>Data acquisition</b>           | Laser power control<br>Exposure time control<br>Sensor gain control   | Number of repetitions control<br>Spectral range control  |
| <b>Data preprocessing options</b> | <b>Spike correction</b><br>(Whittaker-Hayes, moving window)<br><br><b>Spectral smoothing</b><br>(Whittaker, asymmetric least squares, Savitzky-Golay)<br><br><b>Baseline correction</b><br>(Rolling circle, rubberband, least squares, asymmetrically reweighted penalized) | <b>Spectrum normalization</b><br>(Z-score, mean, Mean centre, MinMax, Unit Norm L1, Unit Norm L2)<br><br><b>Intensity normalization*</b><br><br><b>Laser mode-hop correction*</b><br><br><b>Spectral super resolution*</b> |

\*Patented feature

## Data analysis features

|                              |   |  |
|------------------------------|---|--|
| <b>Data Exploration</b>      | Principal component analysis<br>Non-negative matrix factorization (SIMPLISMA-NNLS, MCR-ALS)   |  |
| <b>Qualitative analysis</b>  | Material identification<br>(Pearson correlation, square Euclidean cosine, square first difference Euclidean cosine)<br><br>Creation of spectral libraries |  |
| <b>Quantitative analysis</b> | Principal component regression<br>Segmented HQI option  | Raman peak height/peak area calibration<br>Partial Least Squares (PLS) calibration |
| <b>Classification</b>        | Random Forest, Linear SVM, AdaBoost, Decision Tree, Hoefding Tree, Naïve Bayes, Perceptron, Softmax Regression  |  |

\*Patented feature



## Accessories

- Long working distance probe, f=30 mm
- Middle working distance probe, f=15 mm
- Short working distance probe, f=6.25 mm
- Contact probe
- Disposable silicone heads
- Sample holders
- Axial focusing accessories
- Light protection sample cover
- Power bank
- Belt Holder
- Docking station


## Applications

For researchers, industries, universities and general consumers

- Life Sciences & Health:**  
 Biosciences, pharmaceuticals, skin diagnostics, cosmetics
- Materials & Nanotechnology:**  
 Polymers, nano-materials, semiconductors, surface enhanced Raman scattering (SERS)
- Industrial & Chemical Analysis:**  
 Chemicals, geology, forensics
- Quality Control & Authentication:**  
 Alcohol quality, counterfeit product detection



## Unique Technology



The **miniRaman Spectrometer** is an innovative, miniaturized Raman spectrometer equipped with a patented, integrated **reference channel**. This innovative technology ensures **automatic calibration** by correcting for laser wavelength drift. The system continuously adjusts the Raman shift and intensity of the sample spectrum, delivering fast, accurate, and reliable results. These advanced features make the miniRaman spectrometer an ideal solution for **material identification** and precise quantitative measurements.

## Specification

| Feature versus model     | miniRaman Standard        | miniRaman SERS | miniRaman Dual  | miniRaman Access          | miniRaman Access Full Range |
|--------------------------|---------------------------|----------------|---|---------------------------|-----------------------------|
| Laser wave length        | 785 nm                    |                | 660 nm and 785 nm   | 660 nm                    |                             |
| Power on Sample*         | 5-50 mW                   | 0.5-15 mW      | 5-35 mW (660)<br>5-50 mW (785)                                      | 10-80 mW                  |                             |
| Spectral Range           | 270-3000 cm <sup>-1</sup> |                | 2700-3900 cm <sup>-1</sup> (660)<br>300-2700 cm <sup>-1</sup> (785) | 270-2500 cm <sup>-1</sup> | 270-3900 cm <sup>-1</sup>   |
| Spectral Resolution**    | 8,10,14 cm <sup>-1</sup>  |                | 12,16,20 cm <sup>-1</sup> (660)<br>8,10,14 cm <sup>-1</sup> (785)   | 12,16,20 cm <sup>-1</sup> |                             |
| Signal-to-noise ratio*** | 500:1                     | 100:1          | 800:1 (660)<br>440:1 (785)  | 800:1                     |                             |

\* Actual laser power range can differ  $\pm 2\%$  per device. Please contact us if you need specific laser power range values

\*\* Slit size dependent; slit size can be customized (20,35,50  $\mu$ m slits)

\*\*\* Determined as peak signal-to-noise ratio of polystyrene spectrum at maximal laser power, integration time 0.3s, number of repetitions 10.

## Extra Features

- Small size (105 x 50 x 45 mm)
- Lightweight (360 g)
- One button operation
- Low power consumption (Up to 4 hr)
- Bluetooth & USB connectivity



## About Lightnovo

A spin-off from the Technical University of Denmark, Lightnovo was founded in 2019 by an enthusiastic team united by the goal of revolutionizing the field of Raman spectroscopy through innovative, high-performance solutions. Our mission is to develop and commercialize "Raman for all: democratize the power of high-end Raman spectroscopy for the benefit of mankind".

**We aim to provide premium performance Raman spectrometers and microscopes with the world's smallest form factors without compromising the performance. With this innovation, Lightnovo addresses the need for portable, reliable field instruments at an affordable price.**

## Contact Details

### Lightnovo ApS

Blokken 11, 1.

3460 Birkerød

Denmark (DK)

+45 71 37 04 10

info@lightnovo.com

CVR: 40979603

## Distributor in Ukraine



### LEM Ukraine

Rybalska Str., 2, off. 302

Kyiv, 01011, Ukraine


+380 95 637 78 52

info@lem-ukraine.com.ua

lem-ukraine.com.ua

© 2025 Lightnovo ApS. All rights reserved

Scan to learn more about  
this and other products



Please consult your local sales  
representative for details