

Thermomechanical analysis (TMA) is a technique widely used to determine transition temperatures and expansion/contraction properties of materials, including the coefficient of thermal expansion ( $\alpha$ ). When determining  $\alpha$ , however, there are three factors which need to be considered to obtain reproducible quantitative results. These are TMA sensitivity and baseline drift, sample thickness, and the magnitude of  $\alpha$ .

The TA Instruments TMA has a maximum sensitivity of 15 nm, a displacement resolution of 0.5 nm, and a baseline drift of 1  $\mu\text{m}$  (typically, 0.5  $\mu\text{m}$  between 0 to 500  $^{\circ}\text{C}$ ). Therefore, the recommended measured sample dimensional change over this range for a determinations is >5  $\mu\text{m}$  (roughly 10x the baseline drift). Knowing this, the minimum sample thickness required for testing can be determined from:

$$\alpha = \frac{\Delta L}{\Delta T} * \frac{1}{L_o} \rightarrow L_o = \frac{\Delta L}{\Delta T} * \frac{1}{\alpha} \quad (1)$$

where:

$\Delta L$  = measured dimension change

$L_o$  = original sample thickness

$\Delta T$  = temperature range of measurement

For example, a material with an expected  $\alpha$  of 50  $\mu\text{m}/\text{m}^{\circ}\text{C}$  would yield a recommended minimum sample thickness for testing of:

$$L_o = \left( \frac{5 \mu\text{m}}{500 \text{ }^{\circ}\text{C}} \right) * \left( \frac{1}{50 \mu\text{m}/\text{m} \text{ }^{\circ}\text{C}} \right) = 200 \mu\text{m} \quad (2)$$

This minimum sample thickness will obviously change with a and temperature range. Clearly, it is better to use samples much thicker than this minimum if possible, as this will further improve signal-to-noise ratio.

## ACKNOWLEDGMENTS

This applications note was submitted by Jon Foreman of the Applications Laboratory (US).

For more information or to request a product quote, please visit [www.tainstruments.com/](http://www.tainstruments.com/) to locate your local sales office information.