


## INTRODUCTION

Heat Deflection Temperature (HDT) and Deflection Temperature Under Load (DTUL) are equivalent terms that reflect the temperature at which a material subjected to a three-point bending load deforms to a pre-determined position. The actual force applied to the sample and the amount of deflection required depend upon the sample geometry. ASTM E2092 and a related standard D648, defines a standard set of these parameters. The TMA flexure probe configuration (Figure 1) is utilized in this method.



Flexure Probe

Figure 1.

ASTM E2092 defines DTUL as the temperature at which a precise strain (either 0.25 mm deflection or 0.20% strain as defined by sample dimensions in the procedure\*) occurs under a specific stress (either 455 or 1820 kPa). With the TMA, the load (force) needed to achieve these stresses can be determined using equation (1).

$$F = \frac{2}{3} \frac{Sbd^2}{L} \quad (1)$$

where:

- F = TMA Force (N)
- S = Stress (MPa)  
= 0.455 MPa (66 psi) or 1.82 MPa (264 psi)
- b = Sample width (mm)
- d = Sample thickness (mm)
- L = Sample length = 5.08 mm (as defined by the flexure probe geometry, see Figure 1)

The TMA deflection equivalent to the ASTM specified % strain can be determined using the relationship shown in equation (2).

$$D = \frac{rL^2}{6d} \quad (2)$$

where:

D = TMA dimension change at center span (mm)

r = Sample strain = 0.0020 or 0.20%\*

\*[Sample strain computed based on sample dimensions and deflection change defined by ASTM E2092.

## TMA EXPERIMENTAL GUIDELINES

- Thin samples (1 mm thick or less) are cut to appropriate dimensions (10 mm long x 3 mm wide).
- Samples are placed on the bending fixture and the knife-edged probe is lowered onto the sample with the predetermined force (F).
- The sample is heated at 2 °C/minute until the required deformation is achieved as defined by equation (2).
- The temperature at the required deformation is recorded as the DTUL for the corresponding stress.
- An example of a DTUL experiment is shown in Figure 2 for Polyvinyl Chloride.




Figure 2. Polyvinyl Chloride DTUL - 455 kPa

## REFERENCES

1. ASTM Standard Test Method E2092 & D648.

For more information or to request a product quote, please visit [www.tainstruments.com/](http://www.tainstruments.com/) to locate your local sales office information.